

**Hemin (porcine)**

sterile 500x stock solution in 30% triethanolamine

Cat. No.	Amount
ML-108S	2 ml (for 1 L medium)
ML-108L	10 ml (for 5 L medium)
ML-108XL	20 ml (for 10 L medium)
ML-108XXL	5x 20 ml (for 50 L medium)

**For general laboratory use.** Not intended for human or animal diagnostic or therapeutic uses.

**Shipping:** shipped on gel packs

**Storage Conditions:** store at 4 °C

**Additional Storage Conditions:** store dark

**Shelf Life:** 12 months

**CAS#:** 16009-13-5 (Hemin), 102-71-6 (Triethanolamine)

**EC number:** 240-140-1 (Hemin), 203-049-8 (Triethanolamine)

**Description:**

0,25% (w/v) solution of porcine Hemin in 30% Triethanolamine, sterile. For preparation of growth media for *Leishmania tarentolae*.

**Preparation of *Leishmania* cultivation media:**

Add to 1 l medium:

2.0 ml Hemin stock solution (500x) to a final concentration of 5 µg/ml.