


**1-¹⁵N-rG-CEP**[1-¹⁵N]-N2-Acetyl-5'-O-DMT-2'-O-TBDMS-guanosine-3'-O-CE phosphoramidite

| Cat. No. | Amount |
|------------|--------|
| N-IT-006-S | 10 mg |
| N-IT-006-L | 50 mg |

Applications:

RNA oligonucleotide synthesis (bottle configured for ABI)

Isotopic enrichment: ≥ 98% (Isotopic Enrichment of Precursor)**For general laboratory use.****Shipping:** shipped at ambient temperature**Storage Conditions:** store at -20 °C**Shelf Life:** 12 months after date of delivery**Molecular Formula:** C₄₈H₆₄N₆¹⁵NO₉PSi**Molecular Weight:** 943.13 g/mol**Exact Mass:** 942.42 g/mol**CAS#:** 944138-03-8 (unlabeled)**Purity:** ≥ 95 % (31P-NMR)**Form:** solid**Color:** white of off-white