



EDA-ATP-ATTO-540Q

2'/3'-O-(2-Aminoethyl-carbamoyl)-Adenosine-5'-triphosphate, labeled with ATTO 540Q, Triethylammonium salt

Cat. No.	Amount
NU-808-540Q	40 µl (1 mM)

For general laboratory use.

Shipping: shipped on gel packs

Storage Conditions: store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery

Molecular Formula: C₁₃H₂₂N₇O₁₄P₃ - ATTO 540Q (free acid)

Molecular Weight: 1133.27 g/mol (free acid)



Purity: ≥ 95 % (HPLC)

Form: solution in water


Concentration: 1.0 mM - 1.1 mM

pH: 7.5 ± 0.5

Spectroscopic Properties: λ_{exc} 543 nm, ε 105.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5)



Structural formula of EDA-ATP-ATTO-540Q



absorption spectrum of ATTO 540Q