



## 8-(6-Aminoethyl)-amino-adenosine-3',5'-bisphosphate-DYQ-660

8-(6-Aminoethyl)-amino-adenosine-3',5'-bisphosphate, labeled with DYQ 660, Triethylammonium salt

| Cat. No.    | Amount       |
|-------------|--------------|
| NU-811-660Q | 40 µl (1 mM) |

**For general laboratory use.**

**Shipping:** shipped on gel packs

**Storage Conditions:** store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

**Shelf Life:** 12 months after date of delivery

**Molecular Formula:** C<sub>55</sub>H<sub>71</sub>N<sub>9</sub>O<sub>15</sub>P<sub>2</sub>S (free acid)

**Molecular Weight:** 1192.22 g/mol (free acid)

**Purity:** ≥ 95 % (HPLC)


**Form:** solution in water

**Concentration:** 1.0 mM - 1.1 mM


**pH:** 7.5 ± 0.5

**Spectroscopic Properties:** λ<sub>exc</sub> 660 nm, ε 109.0 L mmol<sup>-1</sup> cm<sup>-1</sup> (Tris-HCl pH 7.5)

λ<sub>exc</sub> 660 nm, ε 140.0 L mmol<sup>-1</sup> cm<sup>-1</sup> (Ethanol)



Structural formula of 8-(6-Aminoethyl)-amino-adenosine-3',5'-bisphosphate-DYQ-660



absorption spectrum of DYQ 660