




5-Propargylamino-CTP-PEG₅-AF647


CTP-PEG₅-AF647

5-(PEG₅-Propargylamino)-cytidine-5'-triphosphate, labeled with AF647, Triethylammonium salt

Cat. No.	Amount
NU-831-PEG5-AF647	10 µl (5 mM)



Structural formula of 5-Propargylamino-CTP-PEG₅-AF647



excitation and emission spectrum of AF647

For general laboratory use.

Shipping: shipped on gel packs

Storage Conditions: store at -20 °C

Short term exposure (up to 1 week cumulative) to ambient temperature possible.

Shelf Life: 12 months after date of delivery

Molecular Formula: C₆₂H₉₀N₇O₃₃P₃S₄ (free acid)

Molecular Weight: 1682.58 g/mol (free acid)

Exact Mass: 1681.37 g/mol (free acid)

Purity: ≥ 95 % (HPLC)

Form: solution in 10 mM Tris-HCl

Color: blue

Concentration: 5.0 mM - 5.5 mM

pH: 7.5 ± 0.5

Spectroscopic Properties: λ_{exc} 649 nm, λ_{em} 671 nm, ε 270.0 L mmol⁻¹ cm⁻¹ (Tris-HCl pH 7.5)

Applications:

RNA Labeling by in vitro Transcription with T7 RNA Polymerase (in-house tested)

Related Products:

HighYield T7 RNA Synthesis Kit, #RNT-101