

**Candida albicans Antigen**

wild strain
native, Solid medium

Cat. No.	Amount
PR-BA117-S	250 µg
PR-BA117-L	1 mg

For general laboratory use.

Shipping: shipped on dry ice

Storage Conditions: store at -80 °C

Additional Storage Conditions: avoid freeze/thaw cycles

Shelf Life: 12 months

Form: liquid (Supplied in 0.1 M glycine buffer pH 9.5)

Culture system: Solid medium

Protein concentration: Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.

Preservation: None

Analysis: ELISA

Preparation: Organisms are washed, disrupted and subjected to an extraction process. The antigen consists of cytoplasmic and cell wall components, presented in alkaline pH glycine buffer.

Inactivation: This product has been inactivated by disruption. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.