

***Helicobacter pylori* Antigen**

native, Solid medium

Cat. No.	Amount
PR-BA118-S	250 µg
PR-BA118-L	1 mg

For general laboratory use.**Shipping:** shipped on dry ice**Storage Conditions:** store at -80 °C**Additional Storage Conditions:** avoid freeze/thaw cycles**Shelf Life:** 12 months**Form:** liquid**Culture system:** Solid medium**Protein concentration:** Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.**Preservation:** None**Analysis:** ELISA**Preparation:** Washed harvested organisms are solubilised and partially purified. Antigen is presented in a 1% non-ionic detergent solution.**Inactivation:** This product has been treated with detergent. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.