



Parainfluenza Virus Type 1 Antigen

Strain Sendai
native, MDCK cell culture

Cat. No.	Amount
PR-BA1261-S	250 µg
PR-BA1261-L	1 mg

For general laboratory use.

Shipping: shipped on dry ice

Storage Conditions: store at -80 °C

Additional Storage Conditions: avoid freeze/thaw cycles

Shelf Life: 12 months

Form: liquid

Culture system: MDCK cell culture

Protein concentration: Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.

Preservation: None

Analysis: ELISA

Preparation: Infected cell culture supernatant is treated to remove cellular material, further purified by size exclusion chromatography and concentrated. The antigen is presented in Tris/NaCl buffer at pH 7.4.

Inactivation: This product has been treated with Beta-propiolactone. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.