



Parainfluenza Virus Type 3 Antigen

Strain C243
native, Vero cell culture

| Cat. No. | Amount |
|-------------|--------|
| PR-BA1263-S | 250 µg |
| PR-BA1263-L | 1 mg |

For general laboratory use.

Shipping: shipped on dry ice

Storage Conditions: store at -80 °C

Additional Storage Conditions: avoid freeze/thaw cycles

Shelf Life: 12 months

Form: liquid

Culture system: Vero cell culture

Protein concentration: Lot specific. Determined by Bradford assay (Biorad) using BSA as standard.

Preservation: None

Analysis: ELISA

Preparation: Infected cell culture supernatant is treated to remove cellular material, further purified by size exclusion chromatography and concentrated. The antigen is presented in Tris/NaCl buffer at pH 7.4.

Inactivation: This product has been treated with Beta-propiolactone. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.