




## C-flat™ 2/1 on 200 copper mesh

Cat. No.	Amount
X-302-CU200	50 pcs.



### Applications:

Support film for biological samples in Cryo-EM techniques

### Description:

20 nm C-flat™ carbon film with 2 µm hole size and 1 µm hole spacing, applied on 200 copper mesh.

**nonsaleable to USA**

**For general laboratory use.**

**Shipping:** shipped at ambient temperature

**Storage Conditions:** store at ambient temperature

**Shelf Life:** n/a